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GB USA

Aesculap® Acculan 3Ti Acculan 3Ti Dermatome GA670

5. Validated reprocessing procedure

Note

Note the validated processing procedure brochure (AVA-V6), available separately, see Validated cleaning and disinfection procedure.

5.1 Single-use products

WARNING

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!
► Do not reprocess the product!

| Designation | Art. no. |
|-----------------|----------|
| Dermatome blade | GB228R |

5.2 Dismantling prior to carrying out the reprocessing procedure

- Remove the batteries from the Acculan 3Ti Dermatome see Removing the battery.
- Open the clamping levers and remove the dermatome blade.
- Disassemble the product immediately after use, as described in the respective instructions for use.

Removing the flap rod of the Acculan 3Ti dermatome

- Loosen nut 2 by turning it clockwise (left-handed thread).
- Screw the nut 2 back to the end of the visible surface.
- Press on nut 2 and push flap rod 18 approx. 4 mm to one side.
- Twist flap rod 18 until it can be removed.
- Remove flap rod 18.
- Slide off the flaps 12 from flap rod 18.

The flap rod 18 has been dismantled, see Fig. 13.

Fig. 13 Flap rod dismantled

Aesculap® Acculan 3Ti Acculan 3Ti Dermatome GA670

5.6 Mechanical cleaning/disinfecting

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

| Phase | Step | T [°C/°F] | t [min] | Water quality | Chemical/Note |
|-------|----------------------|--------------|------------|------------------|---|
| I | Prerinse | <25/77 | 3 | D-W | - |
| II | Cleaning | 55/131 | 10 | FD-W | <ul style="list-style-type: none"> ■ Concentrate, alkaline: <ul style="list-style-type: none"> - pH = 13 - <5 % anionic surfactant ■ 0.5 % working solution <ul style="list-style-type: none"> - pH = 11* |
| III | Intermediate rinse | >10/50 | 1 | FD-W | - |
| IV | Thermal disinfecting | 90/194 | 5 | FD-W | - |
| V | Drying | - | - | - | According to the program for cleaning and disinfection device |

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Exemples

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the equipment and materials applied and personnel involved in the reprocessor control are necessary. Furthermore, the reprocessor should evaluate any deviation from the regarding efficiency and possible adverse consequences.

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